

Week 1: How to load and display a 3D medical image

Materials

MDSC 689.03 - W2017 Students > Image_Datasets

- Hip.zip → MRI image of the hip. DICOM format.
- Thorax.tar.gz → CT image of the thorax. DICOM format.
- head.nii.tar.gz → MRI image of the head. NIfTI format.

*tar.gz files can be uncompressed in linux using ***tar -xvzf *.tar.gz***

*The WinRAR program can also help to uncompress images in Windows and MacOS.

Assignment

- Load images in DICOM format.
- Load images in NIfTI format.
- Display slices of the loaded images.
- **Due on: Next monday at noon. Python files and screen captures must be uploaded to your Dropbox directory.**

Important considerations

- Remember the elements required to display any object in VTK: Reader, mapper, actor, renderer, window, and interactor.
vtkDICOMImageReader, vtkImageMapper, vtkActor2D, vtkRenderer, vtkRenderWindow, and vtkRenderWindowInteractor.
- In order to read DICOM and Nifti images, consider the following readers:
vtkDICOMImageReader and vtkNIFTIImageReader.
- Note the difference between the methods: SetDirectoryName(directory) and SetFileName(file), in each reader.
- Do not forget to set suitable values for the slice you want to display, and the window and level.

```
mapperDicom.SetZSlice(int(zSlice))  
mapperDicom.SetColorWindow(1000)  
mapperDicom.SetColorLevel(0)
```

Extra material

How to navigate through different slices using observers and events.

```
interactor.AddObserver("KeyPressEvent", UserDefinedFunction, 1.0)
```

Some results




Figure 1. Slice number 100 of the Thorax.tar.gz dataset. Window: 1000, level: 0.



Figure 2. Slice number 100 of the head.nii.tar.gz dataset. Window: 1000, level: 500.